

Practical Microbiology - Lecture 7

Biochemical Tests

1) Catalase Test:

This test can be used to detect the enzyme catalase. This enzyme is responsible for protecting bacteria from hydrogen peroxide (H₂O₂) accumulation, which can occur during aerobic metabolism. If hydrogen peroxide accumulates, it becomes toxic to the organism. Catalase breaks H₂O₂ down into water and O₂.



Method:

On the surface of a clean slide mix a few colonies with a drop of hydrogen peroxide, the formation of bubbles indicates a positive test.

It is an essential test to differentiate between the genus *Staphylococcus* and *Streptococcus*, *Staphylococcus* gives a positive result while *Streptococcus* gives a negative result.

2) Oxidase Test:

Oxidase enzymes play an important role in the operation of the electron transport system during aerobic respiration.

Method:

Few drops of freshly prepared oxidase reagent (1% tetramethyl-*para*-phenylene diamine dihydrochloride) are added to a piece of filter paper. A colony from the tested organisms is transferred to the filter paper and rubbed onto the reagent using a wooden stick. A positive reaction is indicated by blue-purple color formation within 1 minute, and a negative reaction by absence of coloration or coloration later than 60 seconds.

This test is useful in the differentiation of *Enterobacteriaceae* from *Pseudomonas* (*Neisseria*, *Vibrio*, and *Pseudomonas*) are oxidase positive.

3- Clumping Factor and Coagulase Test:

Used to differentiate *S. aureus* from other species of the genus *Staphylococcus*.

Clumping Factor: is a surface compound that is responsible for adherence of the organisms to fibrinogen and fibrin. When mixed with plasma, they form clumps. Detection of clumping factor is performed by emulsifying the bacterial colony in a drop of saline on a clean slide to produce a dense uniform suspension, then mixed with a drop of undiluted plasma, a positive reaction is indicated by clumping or aggregation of bacteria within 5-10 seconds.

Coagulase test: coagulase is an enzyme that acts on fibrinogen (soluble plasma protein) and converts it into fibrin (insoluble). This test is performed by tube method by emulsifying a colony of tested bacteria in 1 ml of diluted plasma 1:6 dilution of saline, incubated at 37°C, then examined after (1-4) hours, a positive result is indicated by the clot formation as a result of

conversion the fibrinogen to fibrin, while the negative tubes re-examined after 24 hours.

4) Mannitol Fermentation:

Mannitol is a polysaccharide that is fermented by some bacteria; mannitol semi solid media is prepared in test tubes with adding bromothymole blue as an **indicator**

Principle:

To test the ability of an organism to produce acid end products from mannitol fermentation. Used in identification of *Staphylococcus aureus* which give **positive** test.

Method:

The bacterial samples is inoculated into the media and the media incubated at 37°C for 18 hours, if the sample contain bacteria that fermented mannitol the acid is formed as a results of fermentation, this acid decrease the pH of the media and this will change the color of the indicator from blue to yellow.